# RadiForce®

## RX430 RX430-AR RX430-FC

### Color LCD Monitor

English

Deutsch

Français

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#### Instructions for Use: ■ Setup Manual

**Important:** Please read PRECAUTIONS, this Setup Manual and the User's Manual stored on the CD-ROM carefully to familiarize yourself with safe and effective usage. Please retain this manual for future reference.

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## Compatible Resolutions

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The monitor supports the following resolutions.

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\*<sup>1</sup> Recommended resolution

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\*<sup>2</sup> It is also appropriate for the 10 bit.

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A graphics board in conformance with the VESA standards is required.

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
# RadiForce® RX430 RX430-AR RX430-FC

Color LCD Monitor

## Location of the Caution Statements

Example: RX430

## Symbols on the unit

Symbol	This symbol indicates	
	Main Power Switch:	Press to turn the monitor's main power off.
	Main Power Switch:	Press to turn the monitor's main power on.
	Power button:	Press to turn the monitor's power on or off.
	Alternating current	
	Alerting electrical hazard	
	CAUTION:	Refer to SAFETY SYMBOLS section in PRECAUTIONS.





No part of this manual may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, or otherwise, without the prior written permission of EIZO Corporation. EIZO Corporation is under no obligation to hold any submitted material or information confidential unless prior arrangements are made pursuant to EIZO Corporation's receipt of said information. Although every effort has been made to ensure that this manual provides up-to-date information, please note that EIZO monitor specifications are subject to change without notice.

## Package Contents

Check that all the following items are included in the packaging box. If any items are missing or damaged, contact your local dealer.


### Note

- Please keep the packaging box and materials for future movement or transport of the monitor.
  - Monitor
  - Power cord
  - Digital signal cable (DisplayPort) : PP200
  - Digital signal cable (Dual Link) : DD200-DL
  - EIZO USB cable: MD-C93
  - EIZO LCD Utility Disk (CD-ROM)
  - Setup Manual (this manual)
  - VESA Mounting Screws: M4 x 12, 2 pcs. (RX430/RX430-AR), 6 pcs. (RX430-FC)
  - PRECAUTIONS
  - Cable holder (for RX430-FC only): attachment 1 pc. Clamper 4 pcs. (Large: 3, small: 1) Note: For the installation procedures, refer to page 4.
  - Screws for fixing the cable holder (for RX430-FC only) M4 x 16: 3 pcs., M4 x 8: 2 pcs.
  - Cleaning Kit "ScreenCleaner" (for RX430-FC only).

## Controls and Functions

Example:  
RX430

<b>1. Presence Sensor</b>	Detects the movement of a person in front of the monitor.
<b>2. MODE button</b>	Switches the CAL Switch mode.
<b>3. RETURN button</b>	Cancels the setting/adjustment and exits the Adjustment menu.
<b>4. ▲▼ button</b>	Chooses an adjustment item or increases/decreases adjusted values for advanced adjustments using the Adjustment menu.
<b>5. ENTER button</b>	Displays the Adjustment menu, determines an item on the menu screen, and saves values adjusted.
<b>6. ⏻ button</b>	<ul style="list-style-type: none"> <li>• Turns the power on or off.</li> <li>• Indicates monitor's operation status.</li> </ul> <p>Green: Operating, Orange: Power saving mode, Off: Main power/power off</p>
<b>7. Ambient light sensor</b>	Measures the environmental illumination.
<b>8. Main power switch</b>	Turns the main power on or off.
<b>9. Power connector</b>	Connects the power cord.
<b>10. Input signal connector</b>	DVI-D connector
<b>11. DC5V Output</b>	Supplies electrical power to the external devices.
<b>12. Input signal connector</b>	DisplayPort connector
<b>13. USB port (Up)</b>	Connects the USB cable to use the software that needs USB connection, or to use USB Hub function.
<b>14. Security lock slot</b>	Complies with Kensington's MicroSaver security system.
<b>15. PS/2 connector</b>	Used for maintenance. Do not use for other applications. This may cause the monitor to malfunction.
<b>16. USB port (Down)</b>	Connects a peripheral USB device.
<b>17. Ⓧ mark</b>	Protective earth terminal
<b>18. Stand</b>	Used to adjust the height and angle of the monitor screen.
<b>19. Cable holder</b>	Covers the monitor cables.

\* For instructions for use, refer to the User's Manual on the CD-ROM.

# Connecting Cables

## Attention

- Check that the monitor and the PC are powered off.
- When replacing the current monitor with this monitor, be sure to change the PC settings for resolution and vertical scan frequency to those that are available for this monitor. Refer to the compatible resolution table (back of cover page) before connecting the PC.


## Note

- When connecting multiple PCs to the monitor, refer to the User's Manual on the CD-ROM.

## 1. Connect the signal cables to the input signal connectors and PC.


Check the shape of the connectors, and connect the cables.

After connecting the signal cable, tighten the screws of the connectors to secure the coupling.

## 2. Plug the power cord into a power outlet and the power connector on the monitor.

## 3. Connect the USB cable when using RadiCS LE or ScreenManager Pro for Medical.

## 4. Press to turn on the monitor.

The monitor's power indicator lights up green.

## 5. Turn on the PC.

The screen image appears.


If an image does not appear, refer to "No picture" (page 6) for additional advice.

## Attention

- Turn off the monitor and PC after using them.
- For the maximum power saving, it is recommended that the Power button be turned off. Turning off the main power switch or unplugging the power plug completely shuts off power supply to the monitor.

# Adjusting the Screen Height and Angle

Hold left and right edge of the monitor with both hands, and adjust the screen height, tilt and swivel of the screen to the best condition for working.

## Attention

- Be sure that the cables are correctly connected.
- The stand is not supplied with RX430-FC.

# Attaching the cable holder (RX430-FC)

Attach the cable holder to prevent the cable from falling off.

## 1. Mount the attachment supplied with the monitor.

Mount the attachment using the screws (M4 × 16) while aligning it with the screw holes.


## 2. Hold the clammer on each cable.

Four clammers (large: 3 pcs., small: 1 pc.) are supplied. The small clamper is for the USB cable, and the large clammers are for the power cable and signal cable.

## 3. Fix the clamper.


Mount the clamper using the screws while aligning it with the screw holes.

For the screws to use, refer to the picture shown below.

# Showing Button Guide

Press the front buttons (except  button), the Button Guide appears above the button.


## Note

- The Button Guide will continue to appear while the Adjustment menu or mode menu is showing.
- The Button Guide is displayed differently depending on the selected menu or status.


# Selecting Display Mode


This function allows easy selection of a desired mode according to monitor application.

Mode	Purpose
1-DICOM	Available to display in DICOM mode.
2-Custom	Select one of these user modes to define a desired set of parameter settings.
3-CAL1/4-CAL2/5-CAL3	Displays the screen adjusted by calibration software.
6-Hybrid- $\gamma$	Automatically identifies the display area for the monochrome and color images on the same screen, and displays each set status.
7-sRGB	Suitable for color matching with sRGB compatible peripherals.
8-Text	Suitable for displaying texts for word processing or spreadsheets.

## Procedure

1. Press MODE.  
The mode menu appears at the bottom left of the screen.
2. The mode among the list is highlighted in turn each time MODE is pressed.  
You can switch the mode with  or  while the mode menu is displayed.

# Controlling Monitor Quality

The quality control software "RadiCS LE" (for Windows) is included with this monitor.

To calibrate the monitor and manage the calibration history, install RadiCS LE with reference to the RadiCS LE User's Manual on the CD-ROM. Use the RadiCS SelfQC function and built-in Integrated Front Sensor to perform periodic calibration and Grayscale Check for the monitor independently. For details, refer to the User's Manual on the CD-ROM.

## Attention

- Connect the monitor to the PC with a USB cable when using RadiCS LE.

## Note

- To perform advanced quality control that conforms to medical standards, use the optional monitor quality control software kit "RadiCS UX1".

# No picture

If a problem still remains after applying the suggested remedies, contact your local dealer.


Problem	Possible cause and remedy
<p><b>1. No picture</b></p> <ul style="list-style-type: none"> <li>Power indicator does not light.</li> </ul> <hr/> <ul style="list-style-type: none"> <li>Power indicator is lighting green.</li> </ul> <hr/> <ul style="list-style-type: none"> <li>Power indicator is lighting orange.</li> </ul> <hr/> <ul style="list-style-type: none"> <li>Power indicator flashes orange and green.</li> </ul>	<ul style="list-style-type: none"> <li>Check whether the power cord is connected properly.</li> <li>Turn the main power switch on.</li> <li>Press .</li> <li>Turn off the main power, and then turn it on again a few minutes later.</li> </ul> <hr/> <ul style="list-style-type: none"> <li>Increase &lt;Brightness&gt;, &lt;Contrast&gt;, and/or &lt;Gain&gt; in the Adjustment menu.</li> </ul> <hr/> <ul style="list-style-type: none"> <li>Switch the input signal. For details, refer to the User's Manual on the CD-ROM.</li> <li>Operate the mouse or keyboard.</li> <li>Check whether the PC is turned on.</li> <li>When the Presence Sensor is on, the monitor may have shifted to the power saving mode. Come near to the monitor.</li> </ul> <hr/> <ul style="list-style-type: none"> <li>There is a problem in the device that uses DisplayPort connection. Solve the problem, turn off the monitor, and then turn it on again. Refer to the user's manual of the device that uses DisplayPort connection for details.</li> </ul>
<p><b>2. The message below appears.</b></p> <ul style="list-style-type: none"> <li>This message appears when no signal is input. Example:</li> </ul> <div data-bbox="389 1077 699 1274" data-label="Image"> </div> <ul style="list-style-type: none"> <li>The message shows that the input signal is out of the specified frequency range. (Such signal frequency is displayed in magenta.) Example:</li> </ul> <div data-bbox="389 1422 699 1664" data-label="Image"> </div> <p>fD: Dot clock fH: Horizontal scan frequency fV: Vertical scan frequency</p>	<p>This message appears when the signal is not input correctly even when the monitor functions properly.</p> <ul style="list-style-type: none"> <li>The message shown left may appear, because some PCs do not output the signal soon after power-on.</li> <li>Check whether the PC is turned on.</li> <li>Check whether the signal cable is connected properly.</li> <li>Switch the input signal. For details, refer to the User's Manual on the CD-ROM.</li> </ul> <hr/> <ul style="list-style-type: none"> <li>Check whether the PC is configured to meet the resolution and vertical scan frequency requirements of the monitor ("Compatible Resolutions" in the Setup Manual.).</li> <li>Reboot the PC.</li> <li>Select an appropriate display mode using the graphics board's utility. Refer to the manual of the graphics board for details.</li> </ul>

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