





Protective Screen
FPP5800

Protect your RadiForce 8MP LS560W and LS580W monitor investment with FPP5800 protective screen



- Protects the monitor from front fluid ingress
- Anti-reflective coating for continuous image quality
- Excellent light transmission
- Simple field installation on RadiForce LS560W & LS580W
- Easy cleaning


Accidental damage is not covered by warranty and low quality acrylic sheilds reflect light and degrade image quality. FPP5800 provides reliable protection without sacrificing image quality.



FPP5800 technical specifications

Model Variations (Order No.)	FPP5800 (6GF6200-8AA01-0AA0)
Frame Color	Anthracite Gray
Glass	Laminated Glass
Spectral Transmission	> 98% (at 575 nm wavelength)
Residual Reflection	≤ 1% (at 575 nm wavelength)
Surface Imperfection	DIN 10110-7 5/1x1.0 E 0.63 (100 x 100 mm)
Mounting	Directly mounted on the front of the monitor
Abrasion Resistance	According DIN 58 196-5: 50 strokes (glass only)
Adhesion Test	According DIN 58 196-6-K2 and MIL-C-48497 (glass only)
Salt Spray	According DIN EN ISO 9227 CASS: 24h (glass only)
Net Weight*	13.5 ± 1kg
Operating Temperature	+5 to +45 °C
Storage Temperature	-20 to +60 °C
Storage Humidity	Max. 80 % r. h. non condensing
Operation Humidity	Max. 95 % r. h. non condensing
Operating Pressure	700-1060 hPa (525 -795 mmHg) or up to 3050 m (10,000 ft)
Storage Pressure	500 -1060 hPa (375 -795 mmHg) or up to 5050 m (18,000 ft)
Warranty	2 Years

Dimensions (Unit: mm)



*The additional weight might have impact on the parameters of the ceiling suspension.



EIZO Inc.

5710 Warland Drive, Cypress, California 90630
 Phone: (+1) 562-431-5011, Fax: (+1) 562-431-4811
www.eizo.com/na

All product names are trademarks or registered trademarks of their respective companies. EIZO, RadiForce, FlexScan, ScreenManager, RadiCS and RadiNET are registered trademarks of Eizo Corporation. Specifications are subject to change without notice.